

LogTag® Data Logger

Out of Range Temperature Alarm Guide

Out of Range Temperature Alarms

- The alarm settings should be configured for each device based on the type of storage unit the data logger is used for.
 - *This is why it is important to ensure the device is always connected to the glycol bottle in the correct storage unit.*
 - Reference the device configuration guides for information on how to configure the alarms for the device if they are not already configured.
- NOTE: The alarm does NOT make a sound, it does NOT flash, or do anything other than show the word “ALARM” as shown in the picture below.



- If a temperature goes out of range for more than the configured amount of time, the word “ALARM” will be displayed on the screen until the data is downloaded, regardless of whether the temperature goes back into range.
 - For example: if the refrigerator goes above 8°C from 10 pm – 11 pm and then goes back into range, the alarm will be triggered and the word “ALARM” will be showing on the device when staff arrive in the morning.
- Immediate action must be taken upon discovery of an out of range temperature alarm:
 - Mark the vaccines “DO NOT USE” and continue to store them under proper conditions.
 - *This may require you to move them to your backup location.*
 - Download the data logger information and contact the Immunization Program at 307-777-7952 immediately for further guidance.

Be sure to check out the other LogTag® Resources available at www.immunizewyoming.com!

Contact Kristy Westfall, Compliance Specialist with questions at:
307-777-6519 or Kristy.westfall1@wyo.gov.