Botulinum Toxin

Source

• Clostridium botulinum, a large grampositive, spore-forming, anaerobic bacillus

Characteristics

• Grows anaerobically on Blood Agar and egg yolk plates

Pathogenesis

• Toxin enters nerve terminals and blocks release of acetylcholine, blocking neuro-transmission and resulting in muscle paralysis

Toxicity

- Most lethal of all toxic natural substances
- Groups A, B, E (rarely F) cause illness in humans

Symptoms

- •24-36 h (up to 3 d for wound botulism)
- Progressive skeletal muscle weakness
- Symmetrical descending flaccid paralysis
- Can be confused with stroke, Guillain-Barre syndrome or myasthenia gravis





Gram stain

Lipase on egg yolk plates

Ricin Toxin

Source

• Ricinus communis seeds commonly called castor beans

Characteristics

• Toxin can be disseminated in the form of a liquid, powder or mist

Pathogenesis

- A-chain inactivates ribosomes, interrupting protein synthesis
- B-chain binds to carbohydrate receptors on the cell surface and allows toxin complex to enter cell

Toxicity

- Highly toxic by inhalation, ingestion and injection
- Less toxic by ingestion due to digestive activity and poor absorption
- Low dermal toxicity

Symptoms

- 18-24 h post exposure
- Fever, cough, chest tightness, dyspnea, cyanosis, gastroenteritis and necrosis; death in ~72 h





Ricin plant

Castor beans

Staph Enterotoxin B

Source

• Staphylococcus aureus, a gram-positive cocci

Characteristics

- Appears as grape-like clusters on Gram stain or as small off-white colonies on Blood Agar
- Toxin-producing and non-toxigenic strains of S. aureus will appear morphologically identical

Pathogenesis

 Staphylococcus Enterotoxin B (SEB) is a superantigen. Toxin binds to human class II MHC molecules causing cytokine release and system-wide inflammation

Toxicity

•Toxic by inhalation or ingestion

Symptoms

- 4-10 h post-ingestion, 3-12 h post-inhalation
- •Flu-like symptoms, fever, chills, headache, myalgia
- •Nausea, vomiting, and diarrhea
- Nonproductive cough, chest pain, and dyspnea
- SEB can cause toxic shock syndrome





S. aureus Gram stain

Growth on Blood Agar

Botulinum Toxin

Transmission

- •Aerosol release
- Food contamination
- Wound contamination
- Toxin not transmitted person to person

Clinical Specimens

- Serum: > 5 ml (10 ml preferred)
- Feces and gastric contents: >10 g
- Clinical swab specimens: Place swabs in an anaerobic transport media
- Autopsy: contents from small and large intestines
- Storage: Refrigerate, preferably in plastic containers

Environmental Samples

- •All environmental samples: 100 ml or 2 g
- Food, drinks: Send entire item
- Storage: Refrigerate, preferably in plastic containers

Detection

- Mouse neutralization assay
- •Enzyme-linked immunosorbent assays (ELISA)
- PCR

Ricin Toxin

Transmission

- •Aerosol release
- Food contamination
- Injection
- •Toxin not transmitted person to person

Clinical Specimens

- Serum to test for circulating antibody (>5ml) (Tested at CDC)
- •Urine: >5 ml (10 ml preferred)
- Storage: Room temperature, in plastic containers (do not use glass)
- •Note: Ricin toxin can be denatured by freezing or excess heat

Environmental Samples

- •All environmental samples: 100 ml or 2 g
- •Food, drinks: Send entire item
- Storage: Room temperature, in plastic containers (do not use glass)
- •Note: Ricin toxin can be denatured by freezing or excess heat

Detection

- •Time-resolved fluorescence (TRF)
- PCR
- Ricinine detection by chemical analysis

Staph Enterotoxin B

Transmission

- •Aerosol release
- Food contamination
- •Toxin not transmitted person to person

Clinical Specimens

- 5-10 ml blood in EDTA
- •Urine: > 5 ml (10 ml preferred)
- Respiratory secretions, or nasal swabs
- Bacterial isolates
- Storage: Refrigerate

Environmental Samples

- •All environmental samples: 100 ml or 2 g
- •Food, drinks: Send entire item
- Storage: Refrigerate

Detection

- •Time-resolved fluorescence (TRF)
- •Gel diffusion assay
- Latex agglutination test
- PCR

If any of these toxins are suspected IMMEDIATELY CONTACT the Wyoming Public Health Emergency Response Hot Line at 1-888-996-9104